



## FLEXIBLE ENDOSCOPE LIGHT BUNDLE LIST

Last Update: April 29, 2016

PART NUMBER
11001BN
11001RD1 ES2 (e)
11009 BC 1
11101R ES4
11101RP ES8
11101RP2 ES4
11101RP2 ES8 (e)
11101SK1
11101SK1 ES4
11101SK2 ES4
11272C
11272CU1 (e)
11272CU1 ES2 (e)
11274 BCU1 ES (e)
11274AA ES (e)
11274AAU ES (e)
11275A
11277A
11277AU
11277AU (e)
11278AU1 ES (e)
11278AU1 ES (e2)
11282BN1
11282BN1 ES2
11292ADI ES (e)
11301AB1 ES ( e)
11301BN1
11302BD1 ES2 (e) .016
11302BD2 ES2 (e)
11303BD1
4112-FE (e)
60512-VG
7223-001
7224-001 ES2 (e)
7305-001

PART NUMBER
BF-10
BF-15P
BF-160
BF-160 (e)
BF-1T10
BF-1T130
BF-1T130 (e)
BF-1T160
BF-1T160 (e)
BF-1T180 (e) w-s
BF-1T20
BF-1T200
BF-1T200 (e) w-s
BF-1T20D
BF-1T20D (e)
BF-1T240
BF-1T260
BF-1T30
BF-1T30 (e)
BF-1T40
BF-1T60 (e)
BF-1TQ180 (e) w-s
BF-1TR
BF-20/30 (e) .052
BF-200
BF-20D
BF-240
BF-260
BF-2C30
BF-2T10
BF-30
BF-3C160 ES2 (e) w-s
BF-3C20 (e)
BF-3C30
BF-3C40 (e)

PART NUMBER
A500001A (160")
ACN-1/2
ACN-1/2 ES2
BF-MP160F ES2 (e)
BF-MP60 ES2 (e)
BF-P10/30
BF-P160
BF-P180
BF-P180 ( e ) w-s
BF-P20
BF-P200
BF-P200 (e)
BF-P20D
BF-P20D (e)
BF-P240
BF-P40 (e)
BF-P60 ES2 (e)
BF-PE (e)
BF-Q180 (e) w-s
BF-TE2
BF-XP160F ES (e)
BF-XP40 ES (e)
BF-XT160 (e) w-s
BF-XT20
BF-XT30
BF-XT40
BRO-YL2
BRO-YP2
C-155
C-34
C-8 ES (e)
CF-100HI
CF-100HL
CF-100I
CF-100L
CF-100S
CF-100TI
CF-100TL
CF-10I
CF-10L/20L
CF-10L/CF-20L
CF-130I
CF-130L

PART NUMBER
BF-40
BF-40 (e)
BF-UC160F-OL8 ( e )
CF-140S
CF-1B/W
CF-1T10
CF-1T100I
CF-1T100L
CF-1T10L
CF-1T140L
CF-1T200L
CF-1T20I
CF-1T20L
CF-1TL
CF-200HL
CF-200L
CF-200S
CF-20HL
CF-20I
CF-20S
CF-230I
CF-230L
CF-240AI
CF-240AL
CF-240DL
CF-240I
CF-240L
CF-2T160L
CF-30I
CF-30L
CF-40L
CF-B20
CF-H180AI
CF-H180AL
CF-H190L
CF-H260AL
CF-H260AZL
CF-H260DL
CF-HQ190L
CF-LB3/W
CF-P10L
CF-P10S
CF-P20L

PART NUMBER
CF-130S
CF-140I
CF-140L
CF-Q140L
CF-Q160AI
CF-Q160AL
CF-Q160I
CF-Q160L
CF-Q160S
CF-Q160ZI
CF-Q160ZL
CF-Q165I
CF-Q165L
CF-Q180AI
CF-Q180AL
CF-Q240AL
CF-Q240ZL
CF-Q260AL
CF-Q260DL
CF-V10L
CF-V70
CF-VI
CF-VL
CHF-B160 (e)
CHF-B20
CHF-CB20 ES (e)
CHF-CB30L ES (e)
CHF-CB30S ES (e)
CHF-P10
CHF-P20
CHF-P30
CHF-P30M2
CHF-P60 (e)
CHF-T20 (e)
CYF-2
CYF-2 (e)
CYF-200 (e) w-s
CYF-3
CYF-4
CYF-4 (e)
CYF-5 ES2
CYF-V (e) w-s
CYF-V ES2 (e) w-s

PART NUMBER
CF-P20S
CF-Q140-3200
CF-Q140I
CYF-VHA
DUR-8 ES2
DUR-D
EB-1530T2
EB-1530T3
EB-1570K
EB-1570K ES2
EB-1830
EB-1830T2
EB-1970K
EB-270S
EC-200LR
EC-200M
EC-200MR
EC-201WI
EC-250HL2
EC-250HL5
EC-250LP5
EC-250LS5
EC-250WI
EC-250WI5
EC-250WL
EC-250WL5
EC-250WM
EC-250WM/WI .148.5
EC-250WM5
EC-310LT
EC-3200
EC-3400
EC-3400F
EC-3400L
EC-3430FK
EC-3430L
EC-3430LK
EC-3440F
EC-3440FK
EC-3440L
EC-3470LK
EC-3490FK
EC-3490LI

PART NUMBER
CYF-V2 (e)
CYF-V2 ES2
EC-3800L
EC-3800TL
EC-3801
EC-3801F
EC-3801L
EC-380FKP-FK2P
EC-3830FK
EC-3830FK2
EC-3830L
EC-3830LK
EC-3830MK2
EC-3830TL
EC-3830TLK
EC-3831L
EC-3831LK
EC-3840F
EC-3840FK
EC-3840L
EC-3840T2
EC-3840TL
EC-3870FK
EC-3870FK2
EC-3870LK
EC-3870TLK
EC-3872LK
EC-3890FK
EC-3890LI
EC-3890LK
EC-410LT
EC-450HL2
EC-450HL2 (e)
EC-450HL5
EC-450LP5
EC-450LS5
EC-450WL5
EC-530LS
EC-530WI
EC-530WL
EC-530WM
EC-530WR (e) w-s
EC-590MP (e) w-s

PART NUMBER
EC-3490LK
EC-3490TFI
ECM-3400FE
ECY-1530
ED-200XT (e)
ED-200XT (e)
ED-250XT5 (e)
ED-3410 (e)
ED-3430T
ED-3430TK
ED-3440T (e) w-s
ED-3470TK
ED-3480TK
ED-450XL8 (e)
ED-450XT5 (e)
ED-530XT (e)
EE-1580K
EF-100 (e)
EF-100 ES2 (e)
EF-B114 (e)
EF-B20 (e)
EF-BT28
EF-BT28 (e)
EG-1840
EG-1870K
EG-200FP
EG-200FP/HR
EG-200HR
EG-201FP
EG-2430
EG-250HR
EG-250HR2
EG-250PE5
EG-250WR5
EG-2530
EG-2700
EG-270N5
EG-2730
EG-2730K
EG-2731
EG-2740K
EG-2770K
EG-2790K

<b>PART NUMBER</b>
ECM-3400BS VII
EG-290KP
EG-290P
EG-2930
EG-2930K
EG-2931
EG-2931K
EG-2940
EG-2940K
EG-2970K
EG-2980K
EG-2990I
EG-2990K
EG-300FP
EG-310HR
EG-3300 WL-F
EG-3400
EG-3430
EG-3430K
EG-3440
EG-3630UR
EG-3800T
EG-3830TK
EG-3840T
EG-410CT
EG-410D
EG-450CTS
EG-450HR
EG-450HR2
EG-450WR5
EG-470N5
EG-530WR5 (e)
EN-450P5
EN-450T5
ENF-GP 1/2 Moon
ENF-GP Full Moon
ENF-L3 (e)
ENF-P2 (e)
ENF-P3
ENF-P3 ES3 (e)
ENF-P3/P4 1/2 Moon
ENF-P3/P4 Full Moon
ENF-T3

<b>PART NUMBER</b>
EG-2901 / 2900 Bent(.070)
ENF-V ES2
ENF-V2 ES4
ENF-VH ES4
ENF-VQ
ENF-XP ES
ENT-2000
ENT-3L ES2
EO-270F
ES-3830K
ES-3831
ES-3831K
ES-3870K
EV-G300
EV-G325
FB-10H
FB-10V
FB-10X
FB-15BS
FB-15H 97.5
FB-15P
FB-15RBS
FB-15V
FB-15X
FB-15X ES
FB-18BS ES2
FB-18P
FB-18RBS
FB-18X
FB-19H
FB-19TX
FC-100LR
FC-100MR
FC-300HL
FC-38FX2
FC-38LP
FC-38LX
FCN-15X
FCY-15P2
FCY-15RBS
FEG-201A
FG-100FP
FG-16V

PART NUMBER
FG-1Z
FG-23H
FG-24X
FG-27X
FG-28C
FG-29H
FG-29P
FG-29V
FG-29W
FG-29X
FI-7BS (e)
FI-10BS
FI-10BS (e)
FI-10BS ES
FI-10P
FI-10P/10P2 ES2
FI-10RBS (e)
FI-13BS (e)
FI-13P
FI-16BS ES
FL-100
FL-10P/10P2
FL-16RBS ES
FN-38A
FNL-10P/10P2 ES2
FNL-10RBS
FNL-10S
FNL-13S
FNL-15P2
FNL-15RP3
FNL-15S
FNL-7P2 (e)
FNL-7RP3 (e)
FUR-9P
GF-UC30P ( e )
GF-UE160-AL5 ( e )
GIF-100
GIF-130
GIF-140
GIF-160
GIF-1T100
GIF-1T130
GIF-Q260J

PART NUMBER
GIF-1T140
GIF-1T20
GIF-1T240
GIF-1T30
GIF-1TQ160
GIF-1TV10
GIF-2D
GIF-2T100
GIF-2T160
GIF-2T200
GIF-2T240
GIF-2TQ260M
GIF-E
GIF-FQ260Z
GIF-H180
GIF-H260
GIF-H260Z
GIF-HQ190L
GIF-N180
GIF-N30
GIF-P10
GIF-P140
GIF-P2
GIF-P20
GIF-P230
GIF-P30
GIF-PQ20
GIF-Q
GIF-Q10
GIF-Q140
GIF-Q145
GIF-Q145L
GIF-Q150
GIF-Q160
GIF-Q165
GIF-Q180
GIF-Q20
GIF-Q200
GIF-Q230
GIF-Q240
GIF-Q240X
GIF-Q260
LF-GP ES (e)

<b>PART NUMBER</b>
GIF-Q40
GIF-V
GIF-V10
GIF-V2
GIF-V70
GIF-XP
GIF-XP10
GIF-XP150N
GIF-XP160
GIF-XP180N
GIF-XP20
GIF-XP240
GIF-XP260
GIF-XQ
GIF-XQ10/20
GIF-XQ140
GIF-XQ200
GIF-XQ230
GIF-XQ240
GIF-XQ260
GIF-XQ30
GIF-XQ40
GIF-XT30
GIF-XTQ160
GIF-XV10
HYF-1T (e)
HYF-P (e)
HYF-V ES2 (e)
HYF-XP ES2 (e)
JF-100
JF-130
JF-140
JF-1T20
JF-200
JF-230
JF-240
LF1/LF2
LF-DP/TP (e)
LF-DP/TP ES (e)
LF-DP/TP ES2 (e)
LF-GP (e)
SIF-230
SIF-Q140

<b>PART NUMBER</b>
LF-P (e) Full moon
LF-P ES Full Moon
LF-T (e)
LF-T ES2 (e)
LF-V ES (e)
LF-V ES2 ( e ) w-s
LTF-160
LTF-V2 (e)
LTF-V3 (e)
LTF-VP
MAF-TM ES2 (e)
OGF
OL-1 ES2 (e)
OSF-2
OSF-2U
OSF-3
OSF-V60
PCF-10/20
PCF-100
PCF-130I
PCF-130L
PCF-140I
PCF-140L
PCF-160AI
PCF-160AL
PCF-160AL w/step
PCF-230
PCF-240L
PCF-240S
PCF-H180AL
PCF-H190L
PCF-P240AL
PCF-Q180AI
PCF-Q180AL
PCF-Q260AL
PCF-S
PEF-V (e) w-s
PF-34
SE-1000 (e)
SG22-140
SIF-100

PART NUMBER
SIF-Q180
SIF-Q240
TJF-10/20
TJF-100
TJF-130
TJF-140F
TJF-160
TJF-200
TJF-240
TJF-Q180
TNE-2000BX
U-500
UGI-FP7
URF (e)
URF-P (e)
URF-P2 (e)
URF-P2 ES (e)
URF-P3 ES
URF-P3 ES (e) .018
URF-P3 ES (e) .020
URF-P5 ES
URF-P5 ES (e) .018
URF-P5 ES (e) .020
URF-P6 ES ( e )
URF-V ES2
URF-V ES2 (e)
VB-1530T
VB-1530T2
VB-1830T2
VB-2000
VC-200
VGT-3000
VNL-1170K ES2
VNL-1570TK
VQ-8303A
VR-7142A
VR-7142A
VSF-2900
VSF-3430

PART NUMBER
-------------

**Gulf Fiberoptics Naming Conventions:**

(e) – Indicates that the bundle is made with at least one etched tip.

ES – Indicates that part of the bundle is sheathed with Estane. If this is followed by a number, the number indicates the quantity of legs with Estane sheathing.

Example : ES2 indicates two legs of the bundle have Estane sheathing.

w-s Indicates that one of more etched tips have a stainless steel sleeve, but the etched tip is still exposed.

**CONTACT INFORMATION :**

**Gulf Fiberoptics, Inc.**  
**448 Commerce Blvd.**  
**Oldsmar, Florida 34677 U.S.A.**  
**T: 813-855-6618**  
[www.gulffiberoptics.com](http://www.gulffiberoptics.com)