



T H E G U L F F I B E R O P T I C S G R O U P

July 27, 2011

Dear Valued Customer:

This letter is to inform you that the Gulf Medical Fiberoptics, Inc. fiber optic light cables listed on *Attachment 1* have been validated in the STERRAD® 50, 100S, 200, NX and 100NX Systems. Fiber optic light cables listed in *Attachment 1* that have different fiber bundle sizes and/or different sheathing colors shall also be included.

Testing was performed in accordance with AAMI TIR No.12-2010 guidelines, “Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers.”

Please refer to the STERRAD System User’s Guides for general reprocessing instructions, including proper cleaning, and drying, and packaging information prior to reprocessing any medical device in a STERRAD System. Specific sterilization efficacy questions should be made directly to Advanced Sterilization Products by calling 1-888-STERRAD.

All medical devices should be inspected prior to use to ensure patient safety. Do not use any device if there is a defect or damage noted during the inspection.

It is important to note that Gulf Medical Fiberoptics, Inc. cannot guarantee that any medical device repaired by a third-party company will be compatible in a STERRAD Sterilization System. Please check with the facility for compatibility information prior to reprocessing.

For questions related to the medical device, please contact your local Gulf Medical Fiberoptics representative.

Sincerely,

Craig S. Vogeley
Vice President, Operations

Attachment 1

Device Name	Device p/n	Device Name	Device p/n
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Clip On, 10ft	1102.2101800BK	Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Snap on, 6ft	182.261800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Clip On, 6ft	1102.261800BK	Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Snap on, 8ft	182.281800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Clip On, 8ft	1102.281800BK	Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Olympus, 10ft	192.2101800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: PILLING , 10ft	1112.2101800BK	Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Olympus, 6ft	192.261800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: PILLING , 6ft	1112.261800BK	Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Olympus, 8ft	192.281800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: PILLING , 8ft	1112.281800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF CLIP ON, 10ft	2102.2101800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Clip on, 10ft	1122.2101800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF CLIP ON, 6ft	2102.261800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Clip on, 6ft	1122.261800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF CLIP ON, 8ft	2102.281800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Clip on, 8ft	1122.281800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: PILLING, 10ft	2112.2101800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Snap on, 10ft	162.2101800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: PILLING, 6ft	2112.261800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Snap on, 6ft	162.261800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: PILLING, 8ft	2112.281800BL
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Wolf Snap on, 8ft	162.281800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON CLIP ON, 10ft	2122.2101800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Storz, 10ft	172.2101800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON CLIP ON, 6ft	2122.261800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Storz, 6ft	172.261800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON CLIP ON, 8ft	2122.281800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: Storz, 8ft	172.281800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 10ft	262.2101800BK
Fiber Optic Light Cable - Proximal: ACMI/CIRCON Distal: ACMI/CIRCON Snap on, 10ft	182.2101800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 6ft	262.261800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 8ft	262.281800BK	Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 8ft	262.281800BK

Device Name	Device p/n	Device Name	Device p/n
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 8ft	262.281800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON CLIP ON, 8ft	3122.281800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 10ft	272.2101800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF SNAP ON, 10ft	362.2101800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 6ft	272.261800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF SNAP ON, 6ft	362.261800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: WOLF SNAP ON, 8ft	272.281800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF SNAP ON, 8ft	362.281800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON SNAP ON, 10ft	282.2101800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: STORZ, 10ft	372.2101800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON SNAP ON, 6ft	282.261800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: STORZ, 6ft	372.261800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: ACMI/CIRCON SNAP ON, 8ft	282.281800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: STORZ, 8ft	372.281800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: OLYMPUS, 10ft	292.2101800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON SNAP ON, 10ft	382.2101800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: OLYMPUS, 6ft	292.261800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON SNAP ON, 6ft	382.261800BK
Fiber Optic Light Cable - Proximal: WOLF/DYONICS Distal: OLYMPUS, 8ft	292.281800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON SNAP ON, 8ft	382.281800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF CLIP ON, 10ft	3102.2101800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: OLYMPUS, 10ft	392.2101800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF CLIP ON, 6Ft	3102.261800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: OLYMPUS, 6ft	392.261800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: WOLF CLIP ON, 8Ft	3102.281800BK	Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: OLYMPUS, 8ft	392.281800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: PILLING, 10Ft	3112.2101800BK	Fiber Optic Light Cable - Proximal: STORZ Distal: WOLF CLIP ON, 10ft	4102.2101800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: PILLING, 6Ft	3112.261800BK	Fiber Optic Light Cable - Proximal: STORZ Distal: WOLF CLIP ON, 6ft	4102.261800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: PILLING, 8Ft	3112.281800BK	Fiber Optic Light Cable - Proximal: STORZ Distal: WOLF CLIP ON, 8ft	4102.281800BK

Device Name	Device p/n
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON CLIP ON, 10ft	3122.2101800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON CLIP ON, 6ft	3122.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 8ft	4112.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 10ft	4122.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 6ft	4122.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 8ft	4122.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 10ft	462.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 6ft	462.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 8ft	462.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 10ft	472.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 6ft	472.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 8ft	472.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 10ft	482.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 6ft	482.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 8ft	482.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 10ft	492.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 10ft	573.0101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 6ft	573.061800BK

Device Name	Device p/n
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 10ft	4112.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 6ft	4112.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 6ft	492.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 8ft	492.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 10ft	5102.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 6ft	5102.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 8ft	5102.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 10t	5112.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 6ft	5112.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 8ft	5112.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 10ft	5122.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 6ft	5122.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 8ft	5122.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 10ft	562.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 6ft	562.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 8ft	562.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 8ft	573.081800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON SNAP ON, 10ft	582.2101800BK

Device Name	Device p/n
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON CLIP ON, 10Ft	3122.2101800BK
Fiber Optic Light Cable - Proximal: OLYMPUS/ACMIDistal: ACMI/CIRCON CLIP ON, 6Ft	3122.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 8ft	4112.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 10ft	4122.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 6ft	4122.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON CLIP ON, 8ft	4122.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 10ft	462.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 6ft	462.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: WOLFSNAP ON, 8ft	462.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 10ft	472.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 6ft	472.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: STORZ, 8ft	472.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 10ft	482.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 6ft	482.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: ACMI/CIRCON SNAP ON, 8ft	482.281800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 10ft	492.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 10ft	573.0101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 6ft	573.061800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON SNAP ON, 6ft	582.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: OLYMPUS, 10ft	592.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: OLYMPUS, 6ft	592.261800BK

Device Name	Device p/n
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 10ft	4112.2101800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: PILLING, 6ft	4112.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 6ft	492.261800BK
Fiber Optic Light Cable - Proximal: STORZ Distal: OLYMPUS, 8ft	492.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 10ft	5102.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 6ft	5102.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF CLIP ON, 8ft	5102.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 10t	5112.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 6ft	5112.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: PILLING, 8ft	5112.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 10ft	5122.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 6ft	5122.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON CLIP ON, 8ft	5122.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 10ft	562.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 6ft	562.261800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: WOLF SNAP ON, 8ft	562.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: STORZ, 8ft	573.081800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON SNAP ON, 10ft	582.2101800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: ACMI/CIRCON SNAP ON, 8ft	582.281800BK
Fiber Optic Light Cable - Proximal: PILLING Distal: OLYMPUS, 8Ft	592.281800BK